

-sizing guide- spine

Trend Correx TLSO

Model Num. DCT-5657

PDAC Approved Codes:

OR02- L0456

OR03- L0457

Trend Correx SP TLSO

Model Num. DCT-0464

PDAC Approved Codes:

OR02- L0464

Model Num.	Size	Waist Meas.
DCT-5657; DCT-0464	0	22- 32"
DCT-5657; DCT-0464	1	31- 37"
DCT-5657; DCT-0464	2	36- 40"
DCT-5657; DCT-0464	3	39- 43"
DCT-5657; DCT-0464	4	42- 46"
DCT-5657; DCT-0464	5	45- 49"
DCT-5657; DCT-0464	6	48- 52"
DCT-5657; DCT-0464	7	51- 55"

Notes:

- The above measurement is based on waist circumference
- 10" Extension panel is available for patients needing an excess of 55"