

SIZING GUIDE- KNEE

Trend ROM Hinged Knee

Model Num. DCT-002

Model Num. DC-032

PDAC Approved Codes:

OR02- L1832

OR03- L1833

Model Num.	Size	Meas. (6" above center of patella)
DCT-002-1	XS	Up to 15"
DCT-002-2	S	15- 18"
DCT-002-3	M	18- 21"
DCT-002-4	L	21- 24"
DCT-002-5	XL	24- 27"
DCT-002-6	2XL	27- 30"
DCT-002-7	3XL	30- 35"
DC-032	Universal	Up to 28"