

SIZING GUIDE- KNEE

Trend Pro Fit Hinged Knee

Model Num. DCT-1820

Model Num. DCT-1820U

PDAC Approved Codes:

OR02- L1820

Model Num.	Size	Meas. (6" above center of patella)
DCT-1820-1	XS	Up to 15"
DCT-1820-2	S	15- 18"
DCT-1820-3	M	18- 21"
DCT-1820-4	L	21- 24"
DCT-1820-5	XL	24- 27"
DCT-1820-6	2XL	27- 30"
DCT-1820-7	3XL	30- 35"
DCT-1820U	Universal	Up to 28"