

SIZING GUIDE- KNEE

Trend Pro Fit Hinged Knee

Model Num. DCT-1012

Model Num. DCT-1012U

PDAC Approved Codes:

OR02- L1810

OR03- L1812

Model Num.	Size	Meas. (6" above center of patella)
DCT-1012-1	XS	Up to 15"
DCT-1012-2	S	15- 18"
DCT-1012-3	M	18- 21"
DCT-1012-4	L	21- 24"
DCT-1012-5	XL	24- 27"
DCT-1012-6	2XL	27- 30"
DCT-1012-7	3XL	30- 35"
DCT-1012U	Universal	Up to 28"