

-sizing guide- spine

Trend Quick Form LSO (37)

Model Num. DCT-QF37

PDAC Approved Codes:

OR02- L0637

OR03- L0650

Trend Quick Form (31)

Model Num. DCT-QF31

PDAC Approved Codes:

OR02- L0631

OR03- L0648

Trend Quick Form (27)

Model Num. DCT-QF27

PDAC Approved Codes:

OR02- L0627

OR03- L0642

Model Num.	Size	Waist Meas.
DCT-QF27; QF31; QF37	S	22- 32"
DCT-QF27; QF31; QF37	M	33- 37"
DCT-QF27; QF31; QF37	L	38- 42"
DCT-QF27; QF31; QF37	XL	43- 48"
DCT-QF27; QF31; QF37	2XL	49- 54"
DCT-QF27; QF31; QF37	3XL	55- 60"

Notes:

- The above measurement is based on waist circumference
- 10" Extension panel is available for patients needing an excess of 60"