

Correx TLSO

DCT-5657 The Trend Correx TLSO is a dual-purpose orthosis that supports, stabilizes, and corrects the mid and lower spine (T9-S1). Additionally, the Correx can be used to provide unmatched rotational stability and motion restriction during the post-operative period.

Guided Postural Reeducation

The Trend Correx TLSO safely assists a prolonged low load thoracic stretch to increase spinal flexibility



Highlights • Effectively assist postural correction and reeducation • Provides rotational stability and motion restriction • Low profile design featuring Trend TX2 Breathe Thru Mesh Belt and Trend TX2 Moisture Wicking Lining • Trend TX2 No Grip closure for patients lacking hand strength

Indications • Conditions requiring immobilization of the thoracic and lumbar spine • Spinal Stenosis • Herniated Disc • Degenerative Spinal Fractures • Kyphosis

Patent Pending Dual Pulley System

Telescopes the thoracic core independently from the lumbar brace; this revolutionary mechanism secures a neutral spine alignment

PDAC Approved: L0456/L0457

W trend-med.com P 800.441.2591 E info@trend-med.com

TRENDTM
M E D I C A L